


NC807C Ceiling / Wall Pull Unit

DNU8070100 Rev 6

TYPICAL WIRING CONFIGURATION



Several ceiling / wall pulls may be connected in parallel on one zone.


Calls made from ceiling / wall pulls can be reset at:

- A call point connected in parallel to the ceiling / wall pull
- A reset point connected in parallel to the ceiling / wall pull
- An NC941, NC942 or NC944 call controller
- Any standard indicator panel with its reset button fitted.

Note: Resetting calls at indicator panels is not recommended on nurse call systems as there is a definite requirement for calls to be reset at source.

INSTALLATION GUIDELINES, HINTS AND TIPS

- Always assess the construction of the mounting surface before installation and fix the NC807C to a ceiling or wall using suitable screw fixings.
- Tie the pull cord and its bangles to the ceiling / wall pull as detailed on the pull cord's packaging.
- Keep the cord and bangles clean and tidy by cutting the right hand corner of the bag they are supplied in, threading through 30cm of cord and following steps 1 to 3 below. The cord and bangles can now be left dangling until the bangles are ready to be fitted and the system commissioned.




1




Disassemble the unit and thread the cord through the centre eye.

2



Tie a secure stop knot using the actuator and trim off any excess cord.

3



Reassemble the unit and pull the cord tight.

4

Set the bangles as shown by knotting the cord. Trim off excess.

800-1000mm
above floor

100mm
above floor



Ceiling / wall pulls ARE NOT recommended for applications where there is a risk of ligature (e.g. mental health facilities). In such applications we recommend the use of alternative devices such as tamper-resistant call points. Contact your distributor for details.